

BMP Cytology Funnels for the Shandon CytoSpin® Centrifuge

BMP Cytology Funnels are fully compatible with the Shandon CytoSpin Centrifuge. These high quality Cytology Funnels can be used to deposit a thin layer of cells in a clearly defined area of a microscope slide. The filter card absorbs any excess fluid. Operation is identical to those supplied by Shandon.



BMP Single Cytology Funnel with White Filter Card

- For sample volumes of up to 0.4 mL or less
- Provides consistent and homogeneous cell deposition area of 6 mm
- White Filter Card is pre-attached and is designed for single use only
- Funnel sealing caps included

Part Number	Description	Quantity
BMP-CYTO-S50	Single Cytology Funnel with White Filter Card	50/Pack
BMP-CYTO-S500	Single Cytology Funnel Bulk Pack	500/Case
BMP-CYTO-S50BP	Single Cytology Funnel with Brown Filter Card	50/Pack
BM-FUNNEL-CLIP	Metal Slide Clips	6/Pack



BMP Double Cytology Funnel with White Filter Card

- For sample volumes up to 0.25 mL per slide
- Provides two 6 mm cell deposition areas on each slide
- White Filter Card is pre-attached and is designed for single use only
- Funnel sealing caps included

Part Number	Description	Quantity
BMP-CYTO-DB25	Double Cytology Funnel with White Filter Card	25/Pack
BMP-CYTO-DB500	Double Cytology Funnel Bulk Pack	500/Case
BM-FUNNEL-CLIP	Metal Slide Clips	6/Pack



 **BIOMEDICAL
POLYMERS, INC.**

plastic products for laboratory & medical specimen handling

42 Linus Allain Avenue • Gardner, MA 01440 USA

Toll Free: 800/253-3684 • Tel: 978/632-2555 • Fax: 978/632-2524

Email: sales@biomedicalpolymers.com • www.biomedicalpolymers.com