

# Cytology Funnels with Clip Options

A Product Line of Biomedical Polymers



## BMP Cytology Funnels for the Shandon CytoSpin® Centrifuge:

BMP Cytology Funnels are fully compatible with the Shandon CytoSpin Centrifuge. These high quality Cytology Funnels can be used to deposit a thin layer of cells in a clearly defined area of a microscope slide. The filter card absorbs any excess fluid. Operation is identical to those supplied by Shandon.

- White Filter Card is pre-attached and is designed for single use only
- Funnel sealing caps included



## BMP Single Cytology Funnel with White Filter Card

- For sample volumes of up to 0.4 mL or less
- Provides consistent and homogeneous cell deposition area of 6 mm

Item #	Description	Quantity
BMP-CYTO-S50	Single Cytology Funnel with White Filter Card	50 pack
BMP-CYTO-S500	Single Cytology Funnel Bulk Pack	500 case
BMP-CYTO-S50BP	Single Cytology Funnel with Brown Filter Card	50 pack



## BMP Double Cytology Funnel with White Filter Card

- For sample volumes up to 0.25 mL per slide
- Provides two 6 mm cell deposition areas on each slide

Item #	Description	Quantity
BMP-CYTO-DB25	Double Cytology Funnel with White Filter Card	25 pack
BMP-CYTO-DB500	Double Cytology Funnel Bulk Pack	500 case



## Cytology Funnel Clip Options

- Stainless steel metal clips will not rust or corrode during decontamination.
- Plastic disposable clips save you time.

Item #	Description	Quantity
BMP-FUNNEL-CLIP	Metal Slide Clip	6 pack
BMP-CYTO-PLCLIP	Disposable Cytology Funnel Clip	50 pack

