

## All Tissue Grinders:

- Fully assembled
- Gamma sterilized
- Easy-to-open Tyvek® package
- Includes extra cap
- Instructions For Use
- Includes patient identification label



## BMP Tissue Grinders

### Closed System Prevents Exposure to Sample Aerosols

BMP Tissue Grinders are engineered for safety and ease of use. The closed system design provides the ability to grind the specimen in a sealed container, helping to prevent any risk of contact with the sample. The narrow pestle eliminates the problem of sample overflow common with other tissue grinders. The molded-in abrasive surface does not leave any sediment to obstruct the sample during examination.

Item #	Description	Case Weight	Case Pack
BMP-TG-700	15 mL conical	3 lbs. 1.36 kg	10 per case
BMP-TG-702	35 mL free-standing	4 lbs. 1.81 kg	10 per case
BMP-TG-704	50 mL conical	4 lbs. 1.81 kg	10 per case



Standard Molded-In Abrasive Tip

## BMP Ultra Tissue Grinders

### Vitrified, Silica Cast Pestle Tips Provide Enhanced Grinding Capability for Difficult Samples

The pestle tips of the new BMP Ultra Tissue Grinders have been vitrified using a proprietary silica casting process. The result is a glass-like abrasive surface that provides a more aggressive grinding action to enable efficient processing of difficult and fibrous samples. BMP Ultra Tissue Grinders incorporate the same closed system as our standard grinders.

Item #	Description	Case Weight	Case Pack
BMP-TG-706	15 mL conical	3 lbs. 1.36 kg	10 per case
BMP-TG-708	50 mL conical	4 lbs. 1.81 kg	10 per case



Vitrified Tip

For more information or samples, contact us:

800.253.3684 | bmpmedical.com | sales@bmpmedical.com